

## TeloCol, Acid Soluble Bovine Collagen Solution, 6 mg/mL

Product Number 5026-50ML

TeloCol bovine collagen solution is derived from an acid extraction process yielding a telopeptide-intact collagen. The pro-peptide regions at both ends of the collagen chain, N- and C-telopeptide regions, are maintained.

TeloCol collagen is approximately 95% Type I collagen with the remainder being comprised of Type III collagen. This product is supplied at a concentration of approximately 6 mg/ml with the concentration for each specific lot provided on the product label and Certificate of Analysis that is available with the purchase of each product. TeloCol is sterile filtered and provided in user-friendly packaging for use and storage.

TeloCol is ideal for coating of surfaces, providing preparation of thin layers for culturing cells, or use as a solid gel.

### Specifications

Product Number	5026-50ML
Feature	Acid Soluble Bovine Collagen Solution, Type 1, 6 mg/mL
Size	50 ml

#### CellSystems GmbH

Langer Ring 5 · D-53842 Troisdorf

Phone +49 (0) 2241.25515 - 0 · Fax +49 (0) 2241.25515 - 30  
Internet <http://cellsystems.de> · Email [info@cellsystems.de](mailto:info@cellsystems.de)